



Information for patients attending an MRI Scan with an Abbott FreeStyle Libre Continuous Glucose Monitor

You have been asked to attend for an MRI scan. As you are diabetic and have an Abbott FreeStyle Libre Continuous Glucose Monitor (CGM) there are some risks that you need to be aware of.

Background

Abbott FreeStyle CGM sensors (the **Libre 2, Libre 2 Plus, Libre 3, and Libre 3 Plus**) are safe to scan in the United States (US) under specific MRI scanning conditions. The US sensors are the same as the sensors sold and prescribed in the UK, but the MRI safety labelling is yet to be updated in the UK and Europe.

Given the guidance in the US, NHS Greater Glasgow and Clyde policy is to allow patients with Abbott FreeStyle Libre 2/2+ and 3/3+ sensors to have an MRI Scan following the US conditions for safe MRI scanning.

Risks

You should be aware of the following potential complications with your scan.

- There is a low risk that the FreeStyle Libre sensor could dislodge during your MRI scan.
- There is also a risk that the sensor could stop working due to being exposed to a strong magnetic field.
- MRI scans can cause the thermal heating of implants and there is a potential risk of heating and, or a burn. However, your MRI scan settings will follow the US conditions for safe scanning to reduce this risk as much as possible. We consider the risk of this scenario to be low.

Managing your diabetes

You should be aware that **the FreeStyle Libre sensor will provide inaccurate readings during and up to one hour after your MRI scan**. You should not use the sensor's reported readings to manage your diabetes during this time period. After one hour, the sensor's readings will return to normal, and you can use the sensor readings. If you don't have another method for monitoring your blood glucose levels and this is problematic for you, we can reschedule your MRI scan.

If you use an insulin or infusion pump, you **must** remove this before entering the MRI scanner room.

Further questions

If you have any questions or concerns, please discuss this with your referring clinician or a radiographer.

FreeStyle Libre 2/2+

FreeStyle Libre 3/3+



If you are happy to proceed with your MRI scan, please fill in the information below:

First name(s): _____

Surname: _____

Signature: _____

Date: _____

Author	Owner	Revision	Active Date	Review date	Page
EB	JF	1	23/04/2025	23/04/2028	1 of 1
This document is uncontrolled when printed. Check Revision BEFORE use!					